

Date _____ Deliver to Clinic Date _____

Doctor _____ License Number _____ Phone Number _____

Address _____

Patient Name _____ Gender _____ Age _____

Tooth Number _____ Shade _____ Stump Shade _____

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
32	31	30	29	28	27	26	25	24	23	22	21	20	19	18	17

Items Included with Case

- Master Impression
- Model of Provisional
- Bite Registration
- Analog
- Other _____
- Opposing Impression or Model
- Articulator
- Photos
- Impression Coping
- Diagnostic Wax-Up
- Facebow Transfer Jig
- Photo Card
- Abutment
- Pre-Op Model
- Stick Bite
- Flash Drive
- Final Screw

Type of Restoration Desired

- Porcelain Fused to Metal (PFM)
 - High Noble Noble Chrome Cobalt
- e.Max - Monolithic
- e.Max press over Zirconia
- Prettau (Monolithic Zirconia / MZ)
- Porcelain Butt Margin
- Other _____
- Full Gold Crown (FGC)
 - High Noble Noble
- e.Max - Layered
- Porcelain Fused to Zirconia
- PMMA Temporary (polymethyl methacrylate)
- Diagnostic Wax up

Case Notes and Goals

DOCTOR'S SIGNATURE

Implant Restorations

Screw Retained

- PFM Layered Zirconia e.Max press over Zirconia Implant Temporary

Cement Retained

- Custom Titanium Abutment Custom Zirconia Abutment
 Titanium Hue Gold Hue Zirconia Connection * Metal Connection **
 Custom Cast Gold Abutment Other _____
 Yellow High Noble Yellow Noble

Milling Center Preferred

- Group IV Nobel Biocare Atlantis Straumann 3i No Preference


If Inadequate Clearance:

- Provide Reduction Coping Reduce Prep without Reduction Coping Reduce Opposing Metal/Zirconia Island Please Call Doctor


Return Case for:

- Die Trim Doctor Evaluation after Mounting Framework Try-in Bisque Bake

FRAMEWORK DESIGN



ADDITIONAL SHADE NOTES



* Implants with zirconia connections will be manufactured by Nobel Biocare, Straumann, or Atlantis

**Implants with metal connections on zirconia abutments can be milled by Group IV or Nobel Biocare.

Metal connections are cemented to the zirconia abutment in lab.

Pricing varies by manufacturer, most companies can produce platforms on their competitors, call us or consult your rep for more information. Group IV warranties zirconia custom abutments for breakage or defects in workmanship for 5 years. No warranty is provided for the implant fixture itself. Some implant companies will negate their warranty on the implant if the abutment used is manufactured by a competitor, partner, or third party. Group IV assumes no responsibility for the implant fixture. Please consult your warranty information provided by the implant manufacturer.